




Cryoprobes of 185 mm length are intended for penetrating and freezing soft tissue and can be used in laparoscopic procedures, while shorter cryoprobes are typically used for percutaneous procedures under reduced mechanical stress. It is preferred to use the adapted lengths according to the superficiality of the lesions. Select the cryoprobe based on ice ball size you want to obtain the cooling zone location as defined in the table above.



**Figure 2: For illustration only Color Tag labels for cryoprobe and the matching introducer**



The cryoprobe serial number (S/N) appears on the screen/cryoprobe package/cryoprobe plastic grip.



**Figure 3: The Cryoprobe**

**Cryoprobe connection**




**Cryoprobes are single-use devices supplied in single use sterile packaging.**

**Reuse of single-use devices affects their performance and can cause cross-contamination.**

Connect the cryoprobe to the cryohandle, while maintaining sterility of the cryoprobe:

- 1- Remove the plug that covers the cryoprobe connection point.
- 2- Insert the cryoprobe into the insertion point in the handle as shown on the system screen and screw it until a “NEXT” button appears on screen then confirm screwing by an additional slight rotation to confirm that cryoprobe connection is secured.
- 3- Remove the cryoprobe tip protector.

When done, press **Next** on the screen.



**Figure 4: Cryoprobe connection screen**

**Cryoprobe operation**

Perform a functional **pre-test** to ensure system efficacy and safety as explained in the User Manual.

If a functional problem occurs or there is any unusual appearance (such as frost on the plastic cover near the cryohandle, bubbles or any unusual appearance), press **Cancel** on the system screen and follow system instructions until you are required to safely remove the cryoprobe from the cryohandle. Before operating the

cryoablation system, make sure you have completed all pre-operational stages.



**Do not start the Freeze step before the cryoprobe cool zone center is aligned with the center of the target tissue.**

Before activating the freeze step, insert the cryoprobe into the target tissue under imaging guidance, and follow these steps:

1. Confirm longest dimension of the target tissue
2. Plan the trajectory of the cryoprobe prior to placement. When clinically safe, the center of the cool zone shall be along the longest dimension of the target tissue.
3. The black safety mark should be completely in the tissue.
4. In percutaneous approach, perform a 3 mm skin incision (for example using #11 blade) before the cryoprobe is inserted.
5. Position the tip of the cryoprobe at the distal end of the long axis of the target tissue, when clinically safe
6. Center of the cool zone, should be in the center of the tissue to be ablated (see figure 2 for the distance of the cool zone from the tip of the cryoprobe that is visible under Ultrasound imaging).
7. Maintain sterility and patient safety.

Be aware of the markings on the cryoprobe: the **wide black mark** closest to the tip is the **safety mark**. In percutaneous procedures **it must be completely inside the tissue** to avoid skin burns. The rest of the marks indicate depth of cryoprobe insertion: each mark equals one centimeter with distinctive markings at 5 and 10 cm (50 & 100 mm).



**Figure 5: Illustration of FAP7200000 cryoprobe markings**

Once you have verified that the cool zone center of the cryoprobe is located in the center of the target tissue, you may begin freezing.



**Portions of the cryoprobe other than the freeze zone, including the plastic cover that is located near the cryohandle, may become cold and cause tissue damage. If unwanted freezing occurs, immediately stop the freeze step. In case of frost on shaft, start Extraction if possible. If not, wait for passive Thaw. In both cases, use skin protection techniques.**

**To prevent injury, ice ball growth must be monitored under imaging guidance.**

During **Thaw**, the ice ball melts partially or totally depending on the thaw time and the tissue properties.

Keep the cryoprobe location steady in the target tissue during all of the thaw period. Control the process under Ultrasound or other imaging system.

The **Extraction** step occurs at the end of every treatment. Its purpose is to allow the cryoprobe’s removal from the target tissue in the fastest and safest way.

At the end of the extraction step, a message will be displayed on system screen. Wait for the message, then gently remove the cryoprobe from the target tissue.

If the cryoprobe cannot easily be extracted from the tissue, press **Extraction** on screen to initiate another Extraction step.

**In case the Extraction process isn’t available, wait for passive Thaw.**



After the Extraction step, before extracting the cryoprobe, make sure that the cryoprobe can be easily removed from the tissue.

Do not force removal of the cryoprobe from the tissue as it might increase the risk of tissue damage. Continue the Extraction step or wait for passive thaw until the cryoprobe can be withdrawn easily.

#### Cryoprobe disengagement

After removing the cryoprobe from the target tissue, and only if system screen displays a message that it is safe to disengage the cryoprobe, detach the cryoprobe from the cryohandle as follows:

1. Unscrew the used cryoprobe from the cryohandle and dispose of it appropriately.
2. Remove the single-use sterile cover from the cryohandle.
3. Close the cryohandle with the covering plug.

Following each cryoablation procedure, discard the single use devices (single-use cryoprobe, single-use temperature sensor, cryohandle, flexible hose and touch screen covers and sleeves).

All single use devices are considered to be medical waste and must be disposed of in accordance with medical waste laws and hospital standards. Sharp objects such as the cryoprobe, introducer and temperature sensor must be disposed of in a sharps waste container.

#### Contraindications

Severe infection, uncorrectable coagulopathy, hemodynamic or respiratory instability.

#### Predictable Adverse Events

**General surgery/ Minimal Invasive Procedures - Mild/moderate adverse events:** infection, bleeding, pain, fever, thermal injury, injury to adjacent organs, pneumonia, fall, internal adhesions, incomplete treatment changes in the laboratory parameters-elevation in aspartate aminotransferase and/ or alanine aminotransferase level (this reflects hepatocellular damage), minimal self-limited serum bilirubin level elevation, hemorrhage, hematoma, myoglobinemia, pleural effusion, hemothorax, pneumothorax, thrombosis, diarrhea, nausea, deep vein thrombosis (DVT), transient ischemic attack, hypertension, hypothermia, treatment site reaction, local neuropathy, frostbite, skin burn, vagal reaction, vomiting, needle seeding, user accidental injury. **Severe adverse events:** infection, thermal injury, procedure done on the wrong patient/ part of the body, retention of foreign object after surgery, pulmonary emboli, congestive heart failure, stroke, fall, severe changes in the laboratory parameters-elevation in aspartate aminotransferase and/ or alanine aminotransferase (this reflects hepatocellular damage), life threatening serum bilirubin level elevation, perirenal fluid collection, hemorrhage, cryoshock (hypotension), respiratory compromise, multi organ failure, disseminated intravascular coagulation (DIC), abscess, pleural effusion, hemothorax, pneumothorax, thrombosis of the portal vein branches, allergic/anaphylactoid reaction, angina/coronary ischemia, myocardial infarction, arrhythmia, atelectasis, adjacent organ injury.

**Anesthesia Mild/moderate adverse events:** hypoxia, hypotension -low blood pressure, cough, residual neuromuscular block in recovery, intubation complications, slow to regain consciousness, anesthetic turn off too early.

**Anesthesia Severe adverse events:** hypoxia, hypotension – low blood pressure, residual neuromuscular block in recovery.

**Cryoablation procedure in general - Mild/moderate adverse events:** hematoma, hemorrhage, infection, pain, fever, thermal injury, thrombocytopenia, coagulation dysfunction. **Severe adverse events:** haemorrhage, infection, tumor seeding, thermal

injury, thrombocytopenia, coagulation dysfunction, parenchymal or cryoablated organ injury, incorrect interpretation of post-cryo changes, intra vessels / intra bone gas, emboli, myocardial infarction.

**Percutaneous ablation procedure - Mild/moderate adverse events:** percutaneous hematoma and bleeding infection, adjacent organ injury, CT related adverse effects: radiation, elevated creatinine and renal function injury, reaction. **Severe adverse events:** percutaneous bleeding, infection, adjacent organ injury, **CT related adverse effects:** radiation, acute renal injury, anaphylactic shock in reaction to contrast agent admission, gas emboli.

**Laparoscopic procedure - Mild/moderate adverse events:** vascular and visceral injury, general anesthesia related AE, pneumoperitoneum (that can cause hemodynamic alterations), post procedural abdominal adhesions, abdominal wall hematoma, wound infection and fascial injury. **Severe adverse events:** vascular and visceral injury, general anesthesia related AE, port site metastasis, pneumoperitoneum (that can cause hemodynamic alterations), post procedural abdominal adhesions, cryoprobe and trocar insertion include injuries to major retroperitoneal vessels and to bowel, abdominal wall hematoma, fascial dehiscence and herniation, umbilical hernia, and umbilical wound infection.

**Breast - Fibroadenoma and breast cancer: Mild/ moderate adverse events:** bleeding from breast puncture site, local breast hematoma, local breast infection and thermal injury to the breast skin, breast skin bruising, breast swelling, ecchymosis, edema, fat necrosis. **Severe adverse events:** thermal injury.

**Urology - Renal and Prostate: Mild/moderate adverse events:** adjacent organ injury, allergic reaction, bleeding, DVT, cystitis, Hematuria, hypotension, Idiosyncratic reaction, ileus, infection, pleural effusion, pneumothorax, probe site paresthesia, mild renal failure, renal infarct, voiding dysfunction, urinary tract obstruction, hematomas, ejaculatory dysfunction, erectile dysfunction, penile paresthesia, pelvic pain, perineal pain. **Severe adverse events:** abscess, adjacent organ injury, renal artery/renal vein injury, anaphylactic reaction, coronary ischemia, myocardial infarction, bleeding, death, DVT, hypotension, pulmonary embolism, severe renal failure, stroke, Idiosyncratic reaction, ileus, infection lumbar radiculopathy, pelvic vein thrombosis, pleural effusion, pneumothorax, renal hemorrhage, severe urinary tract obstruction, renal vein thrombosis, rectourethral fistula, scrotal edema, urethral stricture, renal infarct, sepsis, urinary tract leak, urethral sloughing, urethral stricture, urinary fistula, urinary frequency/urgency, urinary incontinence, urinary retention, bladder neck contracture, stroke, need for transfusion due to hemorrhage, tumor seeding, gastro intestinal tract injury.

**Thoracic surgery - Lung and cardiac arrhythmia: Mild/moderate adverse events:** frostbite, respiratory failure, pneumothorax, Hemothorax, pleural effusion, pneumonia, empyema, hemoptysis, lung collapse, thrombosis, phrenic nerve palsy, Pain, fever, cough, back pain, skin injury, pulmonary emboli, Loss of speech (temporary aphasia from laryngeal nerve damage), arm paresis, burn, hemorrhage, pneumonitis. **Severe adverse events:** respiratory failure/arrest, pneumothorax, hemothorax, pleural effusion, pneumonia, empyema, lung collapse, thrombosis, phrenic nerve palsy, subcutaneous emphysema, death due to Acute respiratory distress syndrome, pulmonary emboli, prolonged chest tube drainage, prolonged intubation pulmonary insufficiency/failure, hemorrhage, dyspnea, atelectasis.

**Liver - Mild/moderate adverse events:** changes in the laboratory parameters-elevation in AST (Aspartate Aminotransferase) and/ or ALN (Alanine aminotransferase) level - reflects hepatocellular damage of normal liver parenchyma, minimal self- limited serum bilirubin level elevation, hematoma, hemorrhage, nyoglobinemia, pleural effusion, hemothorax, pneumothorax, thrombosis. **Severe adverse events:** Severe changes in the laboratory parameters-

elevation in AST (aspartate aminotransferase) and/ or ALN (Alanine aminotransferase) level - reflects hepatocellular damage of normal liver parenchyma, life threatening serum bilirubin level elevation, hemorrhage, cryoshock (hypotension, respiratory compromise, multi organ failure, DIC), post procedural abscess- due to conduit ascending bacterial (especially in patients with history of biliary interventions and biliary enteric anastomosis ), pleural effusion, hemothorax, pneumothorax, thrombosis of the portal vein branches.

**Musculoskeletal - Mild/moderate adverse events:** muscular injury, pain, swelling, osteonecrosis, osteomyelitis, chondrolysis, nerve palsy, motor dysfunction, peri-ablational neuropathies. **Severe adverse events:** compartment syndrome osteonecrosis, osteomyelitis, chondrolysis, nerve palsy, motor dysfunction, peri-ablational neuropathies, bowel damage, urinary tract damage, pericardial effusion in (in chest wall ablations), avascular necrosis of femur head, ureteral stricture.

**Gynecology - Mild/ moderate adverse events:** spotting, urinary tract infection, treatment site infection, uterine bleeding, genitourinary perforation, pelvic pain.

**Dermatology - Mild/ moderate adverse events:** skin burn/frostbite, wound complication and wound infection.

**General predictable adverse events when using the introducers:** Trauma to tissue, incomplete treatment, infection, user accidental injury, needle seeding.


**Cryoprobe technical specifications**

The packed cryoprobes shall be stored in a dry, cool, well-ventilated and clean environment without corrosive gas.

In general, IceCure's Cryoprobes are available in various diameters (2.4 mm and 3.4 mm), various ice ball shapes (Spheric, Ellipsoid), various tips (trocar, blunt and pencil) and various lengths (124mm to 185 mm external shaft length) according to the expected application, treated tumor size and surgery approach.

The opening of the cryoprobe pouch should be where the “PEEL HERE” label is positioned.

Figure 6: “Peel Here” label



**Ref number:** FAP 7100000, FAP 7200000, FAP 7400000, FAP 7410000, FAP 7600000, FAP 78000000.

\* Certain configurations are not available in some regions.

**Temperature Range:** -196° C to +40° C  
**Needle diameter:** 2.4 mm (13G) or 3.4 mm (10G)

IceCure Medical, Inc.  
 icecuresupport@icecure-medical.com  
 US/Canada Toll Free 1-888-902-5716

IceCure Medical Ltd.  
 7 HaEshel St., 2<sup>nd</sup> floor,  
 Caesarea 3079504, Israel  
 info@icecure-medical.com  
 Tel: +972-4-623 0333; Fax: +972-4-623 0222



MedNet EC-REP GmbH – Authorized Representative of IceCure Medical Ltd.  
 Borkstraße 10, 48163 Münster, Germany  
 Tel +49 (0) 251 32266-0  
 Fax +49 (0) 251 32266-22